



# Inverted Pendulum Accessory

## A Second Pendulum Option From ECP: Modular System Add-on



### Self-Erecting, Inverted & Noninverted Operation

Lets you control open loop stable and unstable systems and dynamically transition between the two




### Fully Adjustable Dynamic Parameters

Adjustable pendulum weight, rod length, and base inertia are ideal for studying control robustness and supporting multiple student assignments with same equipment



### Easily Adapts to Most ECP Plants





### Two Pendulum Solutions From ECP

ECP now offers two pendulum choices. Our ECP Inverted pendulum (Model 505) is a stand-alone system with unique dynamic characteristics (right half plane zeros and poles). This system has proven itself in 10 years of reliable in-field use and offers one of the lowest cost solutions for a stand-alone experiment commercially available. The Pendulum Accessory (A51) described on this page is based on the classical inverted pendulum model and adds on to other stand-alone ECP systems using their actuators, base feedback, and electronics.

Our inverted pendulum accessory is a cost effective way to enhance your laboratory's capability. It has a wealth of features not found in any other commercially available inverted pendulum apparatus. Precision construction, fully adjustable dynamic parameters, and easy installation make it a valuable addition to any control systems laboratory. As with all ECP systems, complete dynamic models and example controllers are provided along with Matlab® scripts, for easy control modeling and design. The Executive USR® program lets you easily implement control and characterize the system via transient response, frequency response, stability and parameter robustness test features. With provided controllers and plug-and-play installation you will get the system up and running in just minutes and perform interesting experiments the very same day!

### Stimulating Experiments

(Example: High Bandwidth Step Response)



Commanded position (yellow), Base Position (red) and Pendulum Angle (green) test data show classical response characteristics

A high bandwidth stabilizing controller regulates the base position. In the step response, the base makes an initial rapid reverse motion (nonminimum phase controller) [2], then moves in a positive direction [3] and overtakes the rod [4], to decelerate it and come to a final position [5]. The maneuver is complete in 0.8 seconds.